

Laboratory Risk Assessment

<i>Procedure</i>	Diatom Slide Preparation
-------------------------	--------------------------

<i>Level of Risk</i>	Low
-----------------------------	-----

Naphrax slide mountant		
------------------------	--	--

Assessment for C.O.S.H.H.

Procedure	Diatom Slide preparation
------------------	--------------------------

Substance/Procedure	Risk of exposure * L/M/H	HSE Exposure Limits (mg/m ³) (toluene)	Local controls used	Disposal	Emergency procedures
Naphrax Slide mountant --- contains Toluene	L	375	F/C, PPE, DG	F	6, 7